

TECNOSERVICE'21 srl by Tinti Sergio Via Carlo Pisacane, 134 61032 Fano (PU) Italy P.IVA IT0200411413 Tel. +39 0721 805911 FAX +39 0721 809794 e-mail: staff@technochef.it www.technochef.it www.technochef.eu

Services and Technologies for professional catering since 1973



Automatic vacuum thermosealer for bench-top trays:

- made of stainless steel;
- possibility of varying temperature and welding time ;
- seals trays with maximum gastro-norm 1/2 size (mm 325x260);
- availability of molds (not included) with 1 impression for gastro-norm 1/2 trays, with 2 impressions for gastro-norm 1/4 trays and with 4 impressions for gastro-norm 1/8 trays;
- automatic work cycle;
- tank lid with automatic opening ;
- LCD + sensor-VAC digit commands ;
- sensor in the chamber for detecting the percentage of vacuum and gas ;
- inert gas inlet predisposition;
- vacuum pump of 20 meters / cubic / hour ;
- press. max: ± 3 BAR;
- lateral support for transparent film reel;
- manual film dragging;
- ideal for the store and for the laboratory.
- IMPORTANT: to operate, the machine requires an air compressor with a minimum tank of 16/18 liters with an operating pressure of 1.5 atmospheres; Should, in the processing of products such as MEAT ... require for the packaging the use of oxygen-rich gaseous mixtures (percentage greater than 20%) is STRONGLY RECOMMENDED the use of an EXPLOSION-PROOF vacuum pump for safety reasons, available as an option on request.

OPTIONAL - ACCESSORIES :

- 1-mold aluminum mold ;
- 2-mold aluminum mold ;
- 4-mold aluminum mold ;
- Air / gas inlet adjustment pressure gauge ;


CE MARK MADE IN ITALY

	TECHNIC	AL CARD	
	power supply	Monofase	
	Volts	V 230/1	
		5060	
		1,8	
	net weight (Kg)	120	
	breadth (mm)	610	
depth (mm)		600	
	height (mm)	560	
	TECHNIC	AL CARD	
CODE/PICTURES		DESCRIPTION	PRICE/DELIVERY
EM-TRAYGN1/2			Delivery from 4 to 9 days
Skowe (M on Skale) House (M on Skale)	Mod.GN1 / 2 Anticoradal aluminum TRAY800, with 1 impr	TECHNOCHEF - Anticoradal aluminum mold, Mod.GN1 / 2 Anticoradal aluminum mold for Mod.TRAY600 and TRAY800, with 1 impression for Gastro-norm 1/2 containers, 265x320 mm	
EM-TRAYGN1/4			Delivery from 4 to 9 days

EM-TRAYGN1/8



TECHNOCHEF - Anticoradal aluminum mold, Mod.GN1 / 8 Anticoradal aluminum mold for Mod.TRAY600 and TRAY800, with 4 imprints for Gastro-norm trays 1/8, 160x130 mm





Delivery from 4 to 9 days

